

ClearView™ SUB-Q

Subcutaneous Needle Infusion Set Designed Specifically For The Delivery Of Medications To The Subcutaneous Tissue

Enhanced Patient Comfort

Various disk sizes and needles gauges & lengths available to suit the patients needs.

Decreased Risk of Complications

Due to stable anchoring of needle and reduced needle tip movement.

Clear-View disk to observe needle insertion site



The Unique Approach To Accessing The Subcutaneous Tissue Featuring...

Clear, soft insertion/anchoring disk for optimal viewing of insertion site.

Disk lies flat against patient's skin promoting needle stability.

Needle gauges & lengths more suited to patient's body mass.

Transparent dressing included in sterile package.

Improve Compliance through the use of a Dedicated Infusion Set which is ***convenient to use***, is associated with ***less trauma*** and is ***less expensive***.

All with patient satisfaction in mind!

ClearView™ “Sub-Q” Infusion sets

Product List/Ordering Information

Standard Sets

All feature a 90° needle on a clear, flexible anchoring disk for easy visualization of the insertion site. These sets are suitable for most adults with adequate subcutaneous tissue. The sets are offered with a variety of needle and tubing lengths. All infusion and extension sets are individually packaged sterile. Please call if you need a custom set.

Tubing	Needle Lengths		
	6mm	9mm	1.2cm
42'	CVS 43130	CVS 43230	CVS 43330
24"	CVS 33130	CVS 33230	CVS 33330
18"	CVS 23130	CVS 23230	CVS 23330
12"	CVS 63130	CVS 63230	CVS 63330
6"	CVS 03130	CVS 03230	CVS 03330
1" w/o extension	CVS 13130NES*	CVS 13230NES*	CVS 13330NES*
42" Extension Set	CVS 13130	CVS 13230	CVS 13330

Pediatric Sets - The “Skippy” Disk

Designed specifically for the smaller patient with less subcutaneous tissue, they feature a smaller, softer 1/2" clear anchoring disk and shorter 4mm needles. All infusion and extension sets are individually sterile packaged. Please call if you need a custom set.

Tubing	Needle Length	Disk Diameter	Catalog #
42'	4mm	1/2"/Small Disk	CVS 43420
24"	4mm	1/2"/Small Disk	CVS 33420
18"	4mm	1/2"/Small Disk	CVS 23420
12"	4mm	1/2"/Small Disk	CVS 63420
6"	4mm	1/2"/Small Disk	CVS 03420
1" w/o extension	4mm	1/2"/Small Disk	CVS 13420NES*
42" Extension Set	4mm	1/2"/Small Disk	CVS 13420

Notes

All needles are 27 gauge unless specified.

*NES - infusion set does NOT include an extension set.

Transparent dressing are provided in all pouches.

All infusion and extension sets are individually packaged sterile.

Provided in boxes of 8 each / 4 boxes per case.

Special Orders

Any of the above sets can be manufactured with a 24 gauge needle

We can modify any set - please let us know what you need!

7350 N. Ridgeway Avenue Skokie, IL 60076

P 847-674-7075

F 847-674-7066

Norfolk  Medical

www.norfolkmedical.com