

# CLEAR-VIEW “Sub-Q”

## Subcutaneous Needle Infusion Set Designed Specifically For The Delivery Of Medications To The Subcutaneous Tissue

### Enhanced Patient Comfort

*Various disk sizes with needle gauges & length available to suit the patients needs.*

### Decreased Risk of Complications

*Due to stable anchoring of needle and reduced needle tip movement.*



## The Unique Approach To Accessing The Subcutaneous Tissue Featuring...

- ◆ Clear, soft insertion/anchoring disk for optimum viewing of insertion site.
- ◆ Disk lies flat against patient's skin promoting needle stability.
- ◆ Needle gauges & lengths more suited to patient's body mass.
- ◆ Transparent dressing included in sterile package.

Improve Compliance through the use of a Dedicated Infusion Set which is *convenient to use*, is associated with *less trauma* and is *less expensive*.

***All with patient satisfaction in mind!***

**Norfolk  Medical**

7350 North Ridgeway - Skokie, IL - 60076  
Tel: (847) 674-7075 ◆ Fax: (847) 674-7066

# Clear-View “Sub-Q” Infusion Sets

## Product List/Ordering Information

All ClearView “Sub-Q” Sets are packaged 32 sets per case

### *Standard Sets*

All feature a 90° needle on a clear flexible anchoring disk for easy visualization of the insertion site. These sets are suitable for most adults with adequate subcutaneous tissue. The sets are offered with a variety of needle and tubing lengths. All infusion and extension sets are individually sterile packaged. Please call if you need a custom set.

Tubing	Needle Lengths		
	6mm	9mm	1.2cm
42'	CVS 43130	CVS 43230	CVS 43330
24"	CVS 33130	CVS 33230	CVS 33330
18"	CVS 23130	CVS 23230	CVS 23330
12"	CVS 63130	CVS 63230	CVS 63330
6"	CVS 03130	CVS 03230	CVS 03330
1" w/o extension	CVS 13130NES*	CVS 13230NES*	CVS 13330NES*
42" Extension Set	CVS 13130	CVS 13230	CVS 13330

### *Pediatric Sets - The “Skippy” Disk*

Designed specially for the smaller patient with less subcutaneous tissue, they feature a smaller 1/2" clear softer anchoring disk and shorter 4mm needles. All infusion and extension sets are individually sterile packaged. Please call if you need a custom set.

Tubing	Needle Length	Disk Diameter	Catalog #
42'	4mm	1/2" / Small Disk	CVS 43420
24"	4mm	1/2" / Small Disk	CVS 33420
18"	4mm	1/2" / Small Disk	CVS 23420
12"	4mm	1/2" / Small Disk	CVS 63420
6"	4mm	1/2" / Small Disk	CVS 03420
1" w/o extension	4mm	1/2" / Small Disk	CVS 13420NES*
42" Extension Set	4mm	1/2" / Small Disk	CVS 13420

### **NOTES**

All needles are 27 gauge unless specified.

\*NES - infusion set does NOT include an extension set.

Transparent dressing are provided in all pouches.

All infusion and extension sets are individually sterile packaged.

Provided in boxes of 8 each / 4 boxes per case.

### **SPECIAL ORDERS**

Any of the above sets can be manufactured with a 24 gauge needle

We can modify any set - please let us know what you need!