


for **HypoDermoClysis**

An effective alternative to IV hydration therapy.



PRODUCT INFORMATION

For your convenience, all supplies needed for the safe and simple delivery of fluids to the dehydrated, long-term care patient are supplied in a sterile pouch.

Aqua-C Hydration System Product Specification Information

- HDC-2500 High Flow Aqua-C Hydration System** 125cc/hr., Gravity Flow, 36" height
- 2 x (25 gauge x 6mm) needles offset 1.5" on a flexible base with adhesive backing, with a female luer, and 1" tubing
 - 82" 20-drop IV administration set with universal spike and flow regulator
 - prep swab, gauze, latex-free gloves, and transparent dressing
- HDC-2700 Low Flow Aqua-C Hydration System** 80cc/hr., Gravity Flow, 36" height
- 2 x (27 gauge x 6mm) needles offset 1.5" on a flexible base with adhesive backing, with a female luer, and 1" tubing
 - 82" 20-drop IV administration set with universal spike and flow regulator
 - prep swab, gauze, latex-free gloves, and transparent dressing

Aqua-C Hydration Needle Set Product Specification Information

- HDC-2521 High Flow Aqua-C Hydration Needle Set** 125cc/hr., Gravity Flow, 36" height
- 2 x (25 gauge x 6mm) needles offset 1.5" on a flexible base with adhesive backing, with a female luer, and 1" tubing
 - transparent dressing
- HDC-2721 Low Flow Aqua-C Hydration Needle Set** 80cc/hr., Gravity Flow, 36" height
- 2 x (27 gauge x 6mm) needles offset 1.5" on a flexible base with adhesive backing, with a female luer, and 1" tubing
 - transparent dressing